## SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle

Non Pyrogenic/Non Toxic/Latex Free/Rx Only  
Nicht-pyrogen/ungiftig/Latex-frei  
Apyrogène/Non toxique/sans latex  
Apirogénica/no tóxico/libre de látex  
Non pirogenico/Non tossico/Latex Free  
Ikke-pyrogen/Ikke-giftig/Latexfri  
Pyrogeeniton/ei myrkyllinen/lateksiton  
Niet pyrogeen/niet toxisch/latexvrij  
Ikke-pyrogen/ikke giftig/lateksfri  
Apirogena/nietoksyczna/zawartości lateksu  
Apirogénico/Atóxico/Livre de látex

Sol-Millennium Medical, Inc.  
5415 Sugarloaf Parkway, Suite 2203  
Lawrenceville, GA 30043, USA.  
404-973-2200

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## EN Luer Lock Safety Syringe with Exchangeable Needle Instructions for Use

- Peel package as indicated to allow proper removal of SOL-CARE™ Luer Lock Safety Syringe with exchangeable needle and remove needle cap.
- Prepare vial for medication draw up.
- Penetrate center of vial with needle and draw up desired dose of medication. (See Figure 1).
- Prepare injection site and change needle if desired.
- Administer medication per facility protocol.
- To engage the safety feature and deliver the full dose, push the plunger through the zero line. (See Figure 2).
- Pull back on the plunger until the locking ring is engaged. Observe withdrawal of needle into the barrel. (See Figure 3).
- Snap off plunger to permanently disable syringe. (See Figure 4).

### Caution:

- \* This device is sterile unless the package has been damaged or opened.
- \* If unable to activate safety mechanism, discard syringe immediately into an approved sharps container.
- \* This device is single use only.
- \* Maximum needle length cannot exceed 1 1/2" or 38,1mm.

Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

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