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EN Luer Lock Safety Syringe with Exchangeable Needle Instructions for Use

- Peel package as indicated to allow proper removal of SOL-CARE™ Luer Lock Safety Syringe with exchangeable needle and remove needle cap.
- Prepare vial for medication draw up.
- Penetrate center of vial with needle and draw up desired dose of medication. (See Figure 1).
- Prepare injection site and change needle if desired.
- Administer medication per facility protocol.
- To engage the safety feature and deliver the full dose, push the plunger through the zero line. (See Figure 2).
- Pull back on the plunger until the locking ring is engaged. Observe withdrawal of needle into the barrel. (See Figure 3).
- Snap off plunger to permanently disable syringe. (See Figure 4).

Caution:

- * This device is sterile unless the package has been damaged or opened.
- * If unable to activate safety mechanism, discard syringe immediately into an approved sharps container.
- * This device is single use only.
- * Maximum needle length cannot exceed 1 1/2" or 38.1mm.

Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

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